

GRAY MULTIFLEX LONG ARM RETRACTOR

NOVUS | SURGICAL

(OMNI Style Retractor)

Ideal for general, abdominal, pelvic, vascular, orthopedic and urologic procedures



AVAILABLE IN 2 SIZES
Arms can be adjusted
independently

1. NS-GT01 Table Clamp
2. NS-CC01 Coupling Clamp
3. NS-OHA01 Horizontal Arm
4. NS-V02 Long Arm (2 required)
5. NS-EC01 Easy Clamp (6 recommended)
6. NS-V03 End Arm-Hinged (optional)
7. Assortment of retractor blades available

Novus Surgical, Inc.

Office 720-851-0611 | Fax 720-851-0612 | Cell 720-308-6886

10940 S. Parker Rd., Suite 783 | Parker, CO 80134 | www.novus-surgical.com

Advantages of the GRAY Multiflex Long Arm Retractor:

- Double cam action of the Horizontal Arm
- Bevelled Ball provides maximum flexibility
- Arms can be adjusted independently



Recommended Retractor Blades:

- | | |
|-------------|----------------------|
| (4) NS-Y012 | Mayo |
| (4) NS-Y013 | Mayo |
| (1) NS-D042 | Fence/Splanchnic |
| (1) NS-D044 | Fence/Splanchnic |
| (1) NS-K071 | Malleable |
| (1) NS-C080 | Malleable Perforated |
| (1) NS-F061 | Deaver |
| (1) NS-H063 | Deaver |

*Note: Surgeons may choose selection of blades

