


GRAY MULTIFLEX SHORT ARM RETRACTOR

NOVUS | SURGICAL

(OMNI Style Retractor)

Ideal for abdominal, gynaecologic, perineal, and pediatric procedures



AVAILABLE IN 2 SIZES

Arms can be adjusted
independently

1. NS-GT01 Table Clamp
2. NS-CC01 Coupling Clamp
3. NS-OHA01 Horizontal Arm
4. NS-V01i Short Arm (2 required for system)
5. NS-EC03 Easy Clamp (6 recommended) - or -
NS-EC02 Easy Clamps when Mini Blades are used
6. Assortment of retractor blades available

Novus Surgical, Inc.

Office 720-851-0611 | Fax 720-851-0612 | Cell 720-308-6886

10940 S. Parker Rd., Suite 783 | Parker, CO 80134 | www.novus-surgical.com

Advantages of the GRAY Multiflex Short Arm Retractor:

- Double cam action of the Horizontal Arm
- Bevelled Ball provides maximum flexibility
- Arms can be adjusted independently



Recommended Retractor Blades:

- (4) NS-Y009 Mayo Mini
- (4) NS-Y010 Mayo Mini
- (1) NS-D057 Deaver Mini
- (1) NS-D058 Deaver Mini
- (1) NS-K009 Kelly Mini
- (1) NS-C009 Claw Mini
- (1) NS-F009 Finger Mini
- (1) NS-H011 Hoe Mini

*Note: Surgeons may choose any selection of blades

