


NSI CLEAR-LINE LS

Lumbar Retractor System

UPGRADE

YOUR CURRENT LUMBAR RETRACTOR SYSTEM . . .
*with the advantages of
the Clear-Line LS*



Interchangeable Blades:
*with Innovative Blade
Alignment Groove*



**Interchangeable
Frame/Adapters**

NSI CLEAR-LINE LS

BLADE UPGRADE SETS

INTERCHANGEABLE • COLOR-CODED • RADIO-LUCENT



CLEAR-LINE LS COMPLETE – Blade Set

NSI-2335AS	23mm x 35mm	Clear-Line LS; Blade	4	NSI-4665AS	46mm x 65mm	Clear-Line LS; Blade	2
NSI-2345AS	23mm x 45mm	Clear-Line LS; Blade	4	NSI-4675AS	46mm x 75mm	Clear-Line LS; Blade	2
NSI-2355AS	23mm x 55mm	Clear-Line LS; Blade	4	NSI-4685AS	46mm x 85mm	Clear-Line LS; Blade	2
NSI-2365AS	23mm x 65mm	Clear-Line LS; Blade	4	NSI-6055AS	60mm x 55mm	Clear-Line LS; Blade	2
NSI-2375AS	23mm x 75mm	Clear-Line LS; Blade	4	NSI-6065AS	60mm x 65mm	Clear-Line LS; Blade	2
NSI-2385AS	23mm x 85mm	Clear-Line LS; Blade	4	NSI-6075AS	60mm x 75mm	Clear-Line LS; Blade	2
NSI-23100AS	23mm x 100mm	Clear-Line LS; Blade	4	NSI-6085AS	60mm x 85mm	Clear-Line LS; Blade	2
NSI-23120AS	23mm x 120mm	Clear-Line LS; Blade	4	NSI-60100AS	60mm x 100mm	Clear-Line LS; Blade	2
NSI-4655AS	46mm x 55mm	Clear-Line LS; Blade	2	NSI-60120AS	60mm x 120mm	Clear-Line LS; Blade	2

INNOVATIVE



CLEAR-LINE LS BASIC – Blade Set

NSI-2335AS	23mm x 35mm	Clear-Line LS; Blade	4	NSI-2375AS	23mm x 75mm	Clear-Line LS; Blade	4
NSI-2345AS	23mm x 45mm	Clear-Line LS; Blade	4	NSI-2385AS	23mm x 85mm	Clear-Line LS; Blade	4
NSI-2355AS	23mm x 55mm	Clear-Line LS; Blade	4	NSI-23100AS	23mm x 100mm	Clear-Line LS; Blade	4
NSI-2365AS	23mm x 65mm	Clear-Line LS; Blade	4	NSI-23120AS	23mm x 120mm	Clear-Line LS; Blade	4

FRAME UPGRADE SET

CLEAR-LINE LS – Frame Set

NSI-0440	Clear-Line L/LS: Retractor Frame	1
NSI-0437	Clear-Line LS: Side Loading Blade Attachment	6



NSI CLEAR-LINE LS

Proprietary Innovations



Quick-Rail® Retractor Frame *Innovative Curved Rail Design*

- Provides anatomically correct exposure



INNOVATIVE Radiolucent **Clear-Line LS™ Blades**

- Unique blade composition –
Transparent intra-operative scans
- Color-coded



Quick-Lock® Blade Adapters

- Quick / Simple blade attachment
- Works with competitive blades

Optional Items

Quick Release Handles

- Aligns blade for easy
fixation to frame



Clear-Line LS™ **Top Loading** **Blade Attachment**

Novus Surgical, Inc.